A	B	REF	LOT	MD	E	F	G	H	I	J	K	L	M	N	O	P	Q	R	S

REF CNS000004

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Intended use/ Indications: The product is a sterile, single use electrode intended to be attached to the fetal scalp by healthcare professionals during delivery surveillance. The product is intended to be connected to a system for fetal heart-rate monitoring. For use on patients requiring fetal heart rate monitoring during labor. Compatible with: STAN S31 (SBS101003), STAN S41 (SRF618X9)

Patient target group: Fetals during delivery

Intended user: Healthcare professional, typically a midwife or obstetrician, assisting at childbirth.

Contraindication: The FSE should not be applied to fetal face, fontanelles, or genitalia, when placenta previa is present, or in presence of genital infections (e.g. herpes, Group B streptococcus, syphilis, gonorrhoea), Hepatitis, maternal HIV or acquired immune deficiency (AIDS), when mother is a confirmed carrier of hemophilia and fetus is affected or unknown status; or when it is not possible to identify fetal presenting part where application is being considered.

Caution: ⚠ Due to design of fetal electrode, penetration of fetal epidermis may possibly cause trauma, hemorrhage and/or infection. Thus, FSE must be used under conditions of aseptic technique. **AMNIOTIC MEMBRANES MUST BE RUPTURED PRIOR TO ATTACHMENT OF FSE.** Remove FSE from patient before performing any electro-surgical procedures.

Instructions For Use

Application:

- Using aseptic technique, open the package and remove the FSE. Leave Electrode Wires locked in the Retention Notch and check that Protection Tab is in locked position. Straighten the Electrode Wires if necessary.
- Form the Guide Tube gently to the desired anatomical curvature.
- The patient should be in the dorsal lithotomy position. Perform a vaginal examination, to clearly identify the fetal presenting part.
- Insert the FSE at a right angle to the presenting part until the Guide Tube reaches the presenting part. Place the Guide Tube firmly against fetal presenting part.
- To release the Protection Tab, retract the Drive Tube Grip away from the Guide Tube until the Protection Tab clears the Guide Tube (approximately one inch or 2.5 cm).
- Maintain pressure against fetal presenting part with Guide Tube. Advance the Drive Tube until the Spiral Needle reaches the presenting part and with pressure against fetal presenting part, rotate the Drive Tube Grip clockwise until you feel full attachment. **Note:** "full attachment" usually occurs when you feel mild resistance to further rotation and recoil from the Drive Tube Grip. This typically happens after one full rotation.
- Remove the Electrode Wires from the Retention Notch. Then slide Drive Tube of the Electrode Wires. Next slide the Guide Tube of the Electrode Wires and wipe connector clean.
- Insert Connector end securely into Legplate Interface cable.

Removal:

- Pull connector end of FSE out of Legplate.
- Remove Spiral Needle by grasping Electrode Wires as close to fetal presenting part as possible and twist counter-clockwise until free from presenting part. **DO NOT PULL THE SPIRAL NEEDLE FROM THE FETAL SKIN OR PULL WIRES APART.**
- Inspect Spiral Needle tip to make sure it is still attached to Hub. If tip has separated from Hub, remove it from fetal presenting part using aseptic technique.
- Remove Maternal Reference Electrode from Legplate.

All serious incidents occurring in connection with the product should be reported to the manufacturer and the competent authority.

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

A) Consult instructions for use or consult electronic instructions for use B) Catalogue number C) Batch code D) Medical device E) Keep away from sunlight F) Manufacturer G) Date of manufacture H) Use-by date I) Do not use if package is damaged and consult instructions for use J) Humidity limitation 0 – 70% K) Keep dry L) Does not contain or presence of natural rubber latex M) Do not re-use N) Do not resterilize O) Temperature limit 0 - 30°C P) Sterilized using ethylene oxide Q) In accordance with the Medical Device Directive 93/42 / EEC (0402). R) Caution S) Single sterile barrier system

A) Se bruksanvisning eller se elektronisk bruksanvisning B) Katalognummer C) Partikode D) Medisinsk utstyr E) Skal holdes unna sollys F) Produsent G) Produksjonsdato H) Best før-dato I) Skal ikke brukes hvis pakken er skadet, og se bruksanvisning J) Fuktighetsbegrensning 0 – 70 % K) Holdes tørt L) Linneholder ikke naturgummilæte M) Skal ikke gjenbrukes N) Skal ikke steriliseres på nytt O) Temperaturgrense 0 – 30 °C P) Sterilisert med etylenoksid Q) Samsvar med medisinsk utstudsdirektiv 93/42 / EF (0402). R) Forsiktig S) Enkelt sterilt barrieresystem

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Cetro Medical AB, Nitgatan 11, 333 33 Smålandsstenar, Sweden, +46 371 330 30, www.cetromedical.se