

CETRO® Hook

9601 CETRO® Hook

sv en no da fi de fr nl pl et lv lt it el pt es cs
bg ga hr mt ro sk sl hu ar

40-9601-240110

CE
0402

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CETRO® Hook 9601

en

Intended Purpose: The product is a single use, sterile device, used by healthcare personnel for opening the amnion at childbirth, in cases where artificial rupture of the membranes is required.

Intended user: The intended user is a healthcare professional, typically a midwife or obstetrician.

1. Perform a vaginal examination with one hand. Identify the dilation of orifice and the amniotic membrane in front of the presenting fetal part.
2. Hold the cervix aside and gently insert Cetrohook membrane blaster with your other hand.
3. Scratch gently with Cetrohook on the amniotic membrane, regardless of whether the amniotic membrane is in contact with the presenting fetal part or the fetal membrane is bulging.
4. Gently withdraw Cetrohook.

All serious incidents occurring in connection with the product should be reported to the manufacturer and the competent authority.

A) Consult instructions for use **B)** Catalogue number **C)** Batch code **D)** Medical device **E)** Keep away from sunlight **F)** Manufacturer **G)** Date of manufacture **H)** Use-by date **I)** Do not use if package is damaged and consult instructions for use **J)** Humidity limitation 0 – 70% **K)** Keep dry **L)** Does not contain or presence of natural rubber latex **M)** Do not re-use **N)** Do not re-sterilize **O)** Temperature limit 0 - 30°C **P)** Sterilized using ethylene oxide **Q)** In accordance with the Medical Device Directive 93/42 / EEC (0402). **R)** Open here **S)** Single sterile barrier system

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