

Arthroset™

415-2020 Arthroset™ 2,5 mm

415-2025-400 Arthroset™ 4,0 mm

sv en no da fi de fr nl pl et lv lt it el pt es cs
bg ga hr mt ro sk sl hu

40-415-2020/415-2025-400-240503

CE
0402

CETRO MEDICAL

Cetro Medical AB, Nitgatan 11, SE-333 33 Smålandsstenar, Sweden, +46 371 330 30, info@cetromedical.se

415-2020 Arthroset™ 2,5 mm, 415-2025-400 Arthroset™ 4,0 mm

SV

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415-2020 Arthroset™ 2,5 mm, 415-2025-400 Arthroset™ 4,0 mm

en

Intended purpose: The product is a sterile, single use device intended to be used for removal of fluid from an articular joint during arthroscopic procedures, and for the administration of saline solution into the joint. The product is to be used by healthcare professionals.

1. Carefully make sure that the double sterile packing is not damaged. If damage is seen, the cannula set must be discarded and not used in surgery.
2. Open the sterile package by following normal procedure for sterile equipment.
3. Connect the tubing to a relevant syringe used for evacuation of liquid. Make sure the roll clamp is in closed position to avoid spillage.
4. Make a small incision with the tip of a scalpel blade to help easy insertion of the cannula tip.
5. Remove the tip protector from the mandrin's sharp tip (valid for item 415-2020 only).
6. Insert the catheter with the steel mandrin through the incision and place it as required in the operating field.
7. Pull out the steel mandrin from the catheter.
8. Connect the tubing connector to the upper end of the catheter by pressing in carefully into position.
9. For start of fluid outflow, open up the roll clamp to obtain required flow speed.
- 10.If needed, saline solution can be injected into the knee. This is made by attaching the included Luer adapter to the upper part of the catheter. A standard syringe, filled with saline solution, can then be used for filling. Then repeat procedure in pos 8-9.
11. After the procedure, remove the catheter from the knee.

The medical device should be disposed of as per respective hospital policy.

All serious incidents occurring in connection with the product should be reported to the manufacturer and the competent authority.

A)Consult instructions B)Catalogue number C)Batch code D)Medical device E)Keep away from sunlight F)Manufacturer G)Date of manufacture H)Use-by date I)Do not use if package is damaged and consult instructions for use J)Humidity limitation 0 – 70% K)Keep dry L)Does not contain or presence of natural rubber latex M)Do not re-use N)Do not resterilize O)Temperature limit 0 - 30°C P)Sterilized using ethylene oxide Q) In accordance with the Medical Device Directive 93/42/EEC (0402). R)Double sterile barrier system S)Single sterile barrier system

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