
# CETRO® Cordring

9607 CETRO® Cordring 1.0

9607S CETRO® Cordring 2.0


9607

9607S

sv en no da fi de fr nl pl et lv lt it el pt es cs  
bg ga hr mt ro sk sl hu

40-9675-240919

**CE**  
0402

  
**CETRO MEDICAL**

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**Intended purpose /Indications:** The product is a single use, non-invasive, sterile medical device intended to tie off the umbilical cord post-birth.

Application is temporary and made by healthcare professionals.

**Intended user:** The intended user is a healthcare professional, typically a midwife or obstetrician.

**Caution:** See cautions in text below

1. Slip the CETRO® Cordring over the forceps. ⚠ **Do not stretch the CETRO® Cordring more than necessary. This may cause the CETRO® Cordring to break.**
2. Use the cordring forceps to clamp the umbilical cord close to the baby. Ensure that the end tips of the forceps do not project past the diameter of the umbilical cord.
3. Place the cordring forceps next to the other towards the mother. Cut off the umbilical cord between the forceps with a scissor.
4. Securely grasp the handle of the CETRO® Cordring and pull the CETRO® Cordring from the forceps and over the umbilical cord.
5. ⚠ **Do not twist the CETRO® Cordring as this may cause the ring to break.**
6. Cut off the handle of the CETRO® Cordring minimum 1 mm from the ring of the CETRO® Cordring which now is sitting around the umbilical cord. Loosen the forceps.

⚠ **Check the umbilical CETRO® Cordring the first 24 hours after application.**

The medical product must be disposed of according to the respective hospital policy.

All serious incidents that occur in connection with the product must be reported to the manufacturer and competent authority.

A) Consult instructions B) Catalogue number C) Batch code D) Medical device E) Keep away from sunlight F) Manufacturer G) Date of manufacture H) Use-by date I) Do not use if package is damaged and consult instructions for use J) Humidity limitation 0 – 70% K) Keep dry L) Does not contain or presence of natural rubber latex M) Do not re-use N) Do not sterilize O) Temperature limit 0 - 30 °C P) Sterilized using ethylene oxide Q) In accordance with the Medical Device Directive 93/42/EEC (0402). R) Caution S) Latex T) Single sterile barrier system

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