# CETRO® Prolapse Rings


9613 CETRO® Prolapse Ring

9614 CETRO® Prolapse Membrane


9615 CETRO® Prolapse Donut

9616 CETRO® Prolapse Incontinence

9617 CETRO® Prolapse Flexible

sv en no da fi de fr nl pl et lv lt it el pt es cs  
bg ga hr mt ro sk sl hu

40-9605-241011

**CE**  
0402

  
**CETRO MEDICAL**

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en

**Intended use/ Indications:** The product is a sterile, single use, invasive device intended for treatment of uterus prolapse grade I and II, or for treatment of incontinence due to a prolapsed pelvic floor. The product is intended to be inserted into the vagina and gives mechanical support to surrounding tissue.

**Contraindication:** The ring should not be inserted in case of infection in the small pelvis or soft tissue injury in the vagina.

**Caution:** ⚠ Endometriosis can be a reason that prolapse ring should not be used. During pregnancy, the indication is determined on a case-by-case basis.

**Responsibilities:** It is the responsibility of the healthcare professional to inform and provide instruction to the user about treatment, handling, and cleaning. The user shall address all related questions to the responsible healthcare professional.

**BEFORE USE:** -The bladder should be emptied before insertion. -The correct size of the ring is determined by trying it out. -It may require several attempts to try out the right size. -A second visit should suitably take place within 30 days for follow-up symptom reduction and to detect any abrasions in the vagina. -On return visits, the ring should be removed for check-up and the vagina inspected. -The patient should be carefully instructed on the importance of regular check-up by a physician.

**NOTE:** When inserting and removing, the grip should be taken around the ring, not in the holes in the membrane-equipped rings.

- Perform normal gynecological examination of the pelvis and vagina before testing. Use your fingers to appreciate size of ring that may be appropriate to insert.
- Wear dry gloves, lubricate only the front of the ring. Fold the ring according to the recess on the inside of the ring. (Excluding 9615)
- The ring is inserted as shown and opened at the posterior fornix.
- With the index finger, turn the ring a quarter of an hour turns, which prevents the ring from folding and fall out.
- Let the patient stand, sit and bend slightly forward. -A ring that is well placed should not be felt. -The ring should not be too loose nor too tight. One finger should be able to come in between the ring and the vaginal wall. -If there is a large gap between the ring and the vaginal wall it may cause the ring to come out of position or possibly fall out. -When changing to a new ring, do not assume that it will be the same size. -Cetro Medical recommended that you use the same ring for one person for a maximum of 30 days.

The medical product must be disposed of according to the respective hospital policy.

All serious incidents that occur in connection with the product must be reported to the manufacturer and competent authority.

A) Consult instructions B) Catalogue number C) Batch code D) Medical device E) Keep away from sunlight F) Manufacturer G) Date of manufacture H) Use-by date I) Do not use if package is damaged and consult instructions for use J) Humidity limitation 0 – 70% K) Keep dry L) Does not contain or presence of natural rubber latex M) Do not re-use N) Do not sterilize O) Temperature limit 0 - 30°C P) Sterilized using ethylene oxide Q) In accordance with the Medical Device Directive 93/42/EEC (0402). R) Caution S) Single sterile barrier system

9613 CETRO® Prolapse Ring, 9614 CETRO® Prolapse Membrane, 9615 CETRO® Prolapse Donut, 9616 CETRO® Prolapse Incontinence

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